

PRODUCT BROCHURE

KERN PCB-A01

Date: 04.12.2025

ENGLISH



Image not available

Internal rechargeable battery pack, operating time up to 48 h, charging time approx. 8 h, AUTO-OFF function to preserve the rechargeable battery pack

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

Phone : +49 7433 9933-0
Fax : +49 7433 9933-149
Email : info@kern-sohn.com
Web : www.kern-sohn.com



Scan here for more
information

SPECIFICATIONS

Article codes, product groups & Status information

Model article-code	PCB-A01
Model-Serie	PCB
Assortment	KERN
Product Type	Laboratory balances
Product group	Accessories (Printer, etc)
Product usage type	Spare part;Accessories/Services
EAN code (Model)	4045761085657
Predecessor Model	PCB-A06
Customs tariff number	85075000

Technical data - Power supply

Mains adapter/adapter included	CH;EURO;UK;US
Mains adapter	Battery & Power supply

Technical data - Packing & Shipping

Dimensions packaging (W×D×H)	70×100×100 mm
Gross weight	0.40 kg
Net weight	0.25 kg
Shipping method	Parcel service
UN number	3496
Packaging component - by weight - paper	0 g
Packaging component - by weight - cardboard	82.92 g
Packaging component - by weight - plastics	82.92 g

Packaging component - by weight - styrofoam	0 g
Packaging component - by weight - glass	0 g
Packaging component - by weight - wood	0 g
Packaging component - by weight - foam	0 g
Packaging component - by weight - metal	0 g
Delivery time	1 d

FUNCTIONS

Standard

ENGLISH